





Recombinant Mouse RNA binding protein fox-1 homolog 3 (Rbfox3)

Product Code	CSB-EP806847MO
Relevance	Pre-mRNA alternative splicing regulator. Regulates alternative splicing of RBFOX2 to enhance the production of mRNA species that are targeted for nonsense-mediated decay (NMD)
Abbreviation	Recombinant Mouse Rbfox3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8BIF2
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MAQPYPPAQYPPPPQNGIPA EYAPPPHPTQDYS GQTPVPPEHGMTLYTPAQ THPEQPGTEASTQPIAGTQTVPQADEAAQTDNQQLHPSDPTEKQQPKRLHVS NIPFRFRDPDLRQMFGQFGKILDVEIIFNERGSKGFGFVTFETSSDADRAREKL NGTIVEGRKIEVNNATARVMTNKKPGNPYANGWKLNPVVGTVYGPEFYAVTS FPYPTTGTA VAYRGAHLRGRGRAVYNTFRAAPPPPIPTYGAALEQTLVKMPV PWAGLAPCPLPPQQTPEPAYPTSPAFPLSCPFA SRVVYQDGFYGA EIYGGY AAYRYAQPA AATAAA YSDSYGRVYAAADPYHHTIGPTATYSIGTMASLCRGGY SRFTPY
Research Area	Others
Source	E.coli
Target Names	Rbfox3
Protein Names	Fox-1 homolog C (Hexaribonucleotide-binding protein 3) (Fox-3) (Neuronal nuclei antigen) (NeuN antigen) (D11Bwg0517e) (Hrnbp3)
Expression Region	1-374aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	47.6 kDa
Protein Length	Full Length
Image	



We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.