



Recombinant Mouse Protein delta homolog 2 (Dlk2), partial

Product Code	CSB-BP812967MO
Abbreviation	DLK2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8K1E3
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	DDCS SHCDLAHGCC APDGSCRCDP GWEGLHCERC VRMPGCQHGT CHQPWQCICH SGWAGKFC DK DEHICTSQSP CQNGGQCVYD GGGEYHCVCL PGFHGRGCER KAGPCEQAGF PCRNGGQCQD NQG FALNFTC RCLAGFMGAH CEVNVDDCLM RPCANGATCI DGINRFSCLC PEGFAGRFCT INLDDCASRP CQRGARCRDR VHDFDCLCPS GYGGKTCELV LPAPEPASVG TPQMPTSAVV VPATGPAPHS AGAGLLRISV KEVRRQESG LGESS
Source	Baculovirus
Target Names	DLK2
Protein Names	Recommended name: Protein delta homolog 2 Short name= DLK-2 Alternative name(s): Endothelial cell-specific protein S-1 Epidermal growth factor-like protein 9 Short name= EGF-like protein 9
Expression Region	27-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Extracellular domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.