



Recombinant Rabbit Osteopontin (SPP1)

Product Code	CSB-MP022603RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P31097
Product Type	Recombinant Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Purity	>85% (SDS-PAGE)
Sequence	LPVK HADSGSSEEK QLYHKHPDAL ATWLNPDPSQ KQNLLTPQNA MSSEEKDDLK QETLPSKSIE SHDHMDDIDE DEDDDHVDNR DSNESDDADH PDDSHHSDES HQSDESDEVT VYPTEDAATT VFTEVVPTVE TYDGRGDSVA YRLKRSKSKM FHVSNAQYPG ASEEDLSSHV DSEDLDDTPR AIPVAQHLNV PSDWDSQEKD SHDVSQVDDH SVETQSHEQA RQYKREANDN SVEHSHSIDS QESSKVSQES QSREFRSHED KLAIEPKSEE DEEHRQLRVS HELDSTSSEI N
Source	Mammalian cell
Target Names	SPP1
Protein Names	Recommended name: Osteopontin Alternative name(s): Bone sialoprotein 1 OC-1 Secreted phosphoprotein 1 Short name= SPP-1
Expression Region	17-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.