




Recombinant Human Splicing regulatory glutamine/lysine-rich protein 1 (SREK1), partial

Product Code	CSB-EP837871HU
Relevance	Participates in the regulation of alternative splicing by modulating the activity of other splice factors. Inhibits the splicing activity of SFRS1, SFRS2 and SFRS6. Augments the splicing activity of SFRS3
Abbreviation	Recombinant Human SREK1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WXA9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	KERDHISSERTERRERSTSMRKSSNDRDGKEKLEKNSTSLKEKEHNKEPDSSVS KEVDDKDAPRTEENKIQHNGNCQLNEENLSTKTEAV
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	SREK1
Protein Names	Serine/arginine-rich-splicing regulatory protein 86 (SRrp86) (Splicing factor, arginine/serine-rich 12) (Splicing regulatory protein 508) (SRrp508)
Expression Region	421-508aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	17.3 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.