



Recombinant Horse Major allergen Equ c 1

Product Code	CSB-BP839158HO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q95182
Product Type	Recombinant Protein
Immunogen Species	Equus caballus (Horse)
Purity	>85% (SDS-PAGE)
Sequence	QQEEN SDVAIRNFDI SKISGEWYSI FLASDVKEKI EENGSMRVFV DVIRALDNSS LYAEYQTKVN GECTEFPMVF DKTEEDGVYS LNYDGYNVFR ISEFENDEHI ILYLVNFDKD RPFQLFEFYA REPDVSPEIK EEFVKIVQKR GIVKENIIDL TKIDRCFQLR GNGVAQA
Source	Baculovirus
Protein Names	Recommended name: Major allergen Equ c 1 Alternative name(s): Allergen= Equ c 1
Expression Region	16-187
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.