




Recombinant Human Bile acid receptor (NR1H4)


Product Code	CSB-EP839406HU
Relevance	Ligand-activated transcription factor. Receptor for bile acids such as chenodeoxycholic acid, lithocholic acid and deoxycholic acid. Represses the transcription of the cholesterol 7-alpha-hydroxylase gene (CYP7A1) through the induction of NR0B2 or FGF19 expression, via two distinct mechanisms. Activates the intestinal bile acid-binding protein (IBABP). Activates the transcription of bile salt export pump ABCB11 by directly recruiting histone methyltransferase CARM1 to this locus
Abbreviation	Recombinant Human NR1H4 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96R11
Alias	Farnesoid X-activated receptorFarnesol receptor HRR-1Nuclear receptor subfamily 1 group H member 4Retinoid X receptor-interacting protein 14 ;RXR-interacting protein 14
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MGSKMNLIEHSHLPTTDEFSFSENLFGLVTEQVAGPLGQNLEVEPYYSQYSNVQ FPQVQPQISSSSYYSNLGFYPQQPEEWYSPGIYELRRMPAETLYQGETEVAE MPVTKKPRMGASAGRIKGDELVCVCGDRASGYHYNALTCEGCKGFFRRSITK NAVYKCKNGGNCVMDMYMRRKCQECRLRKCKEMGMLAECMYTGLLTEIQCK SKRLRKNVKQHADQTVNEDSEGRDLRQVTSTTKSCREKTELTPDQQTLLHFIM DSYNKQRMPEITNKILKEEFSAEENFLILTEMATNHVQVLVEFTKKLPGFQTL DHEDQIALLKGSAMEFLRSAEIFNKKLPSGHSDDLLEERIRNSGISDEYITPMF SFYKSIGELKMTQEEYALLTAIVILSPDRQYIKDREAVEKLQEPLLDVLQKLCKIH QPENPQHAFACLLGRLTELRTFNHHHAEMLSWRVNDHKFTPLLCEIWDVQ
Research Area	Transcription
Source	E.coli
Target Names	NR1H4
Expression Region	1-476aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	70.7kDa
Protein Length	Full Length of Isoform 3




Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP839406HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) NR1H4.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP839406HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) NR1H4.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.