




Recombinant Human Casein kinase I isoform alpha-like (CSNK1A1L)

Product Code	CSB-EP839816HU
Relevance	Casein kinases are operationally defined by their preferential utilization of acidic proteins such as caseins as substrates. It can phosphorylate a large number of proteins. Participates in Wnt signaling. Phosphorylates CTNNB1 at 'Ser-45'. May phosphorylate PER1 and PER2. May play a role in segregating chromosomes during mitosis.
Abbreviation	Recombinant Human CSNK1A1L protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N752
Alias	CK1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MTNNSGSKAELVVGKYLVRKIGSGSFGDVYLGITTTNGEEVAVKLESQKVK HPQLLYESKLYTILQGGVGIPMHWWYGQEKDNNVLVMDLLGPSLEDLFNFCSR RFTMKTVLMMLADQMISRIEYVHTKNFLHRDIKPDNFLMGTRHCNKLFLIDFGL AKKYRDNRTQRHIPYREDKHLIGTVRYASINAHLGIEQSRRDDMESLG YVFM Y FNRTSLPWQGLKAMTKKQKYEKISEKKMSTPVEVLCKGFPAEFAMYLNYCRG LRFEEVPDYM YLRQLFRILFRTL NHQYDYTFDWTMLKQKAAQQAASSSGQGQ QAQTQTGKQTEKNKNNVKDN
Research Area	others
Source	E.coli
Target Names	CSNK1A1L
Protein Names	Recommended name: Casein kinase I isoform alpha-like Short name= CKI- alpha-like EC= 2.7.11.1Alternative name(s): CK1
Expression Region	1-337aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	66.1kDa
Protein Length	Full Length of BC028723
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.