



Recombinant Bubalus bubalis Follistatin (FST)

Product Code	CSB-BP009024BWX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A5YT95
Product Type	Recombinant Protein
Immunogen Species	Bubalus bubalis (Domestic water buffalo)
Purity	>85% (SDS-PAGE)
Sequence	G NCWLRQAKNG RCQVLYKTEL SKEECCSTGR LSTSWTEEDV NDNTLFKWMI FNGGAPNCIP CKETCENVDC GPGKKCRMNK KNKPRVCAP DCSNITWKGL VCGLDGKTYR NECALLKARC KEQPELQVQY QGKCKKTCRD VFCPGSSTCV VDQTNAYCV TCNRICPEPT SSEQYLCGND GVTYPSACHL RKATCLLGRS IGLAYEGKCI KAKSCDDIQC TGGKKCLWDF KVGRGRCSLC GELCPESKSE EPVCASDNAT YASECAMKEA ACSSGVLLEV KHSGSCNSIS EDTEDEEEDE DQDYSFPISS ILEW
Source	Baculovirus
Target Names	FST
Protein Names	Recommended name: Follistatin Short name= FS Alternative name(s): Activin-binding protein
Expression Region	30-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.