



Recombinant Human Uncharacterized protein C1orf54 (C1orf54)

Product Code	CSB-EP840992HU
Abbreviation	C1orf54
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WWF1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	QEYEDEERLGEDEYYQVVYYYYVTTPSYDDFSADFTIDYSIFESEDRLNRLDKDI TEAIETTISLETARADHPKPVTVKPVTTEPSPDLNDAVSSLRSPIPLLLSCAFVQ VGMYFM
Research Area	Signal Transduction
Source	E.coli
Target Names	C1orf54
Protein Names	Recommended name: Uncharacterized protein C1orf54
Expression Region	17-131aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-B2M-tagged
Mol. Weight	27.2kDa
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.