




# Recombinant Human Phosphoinositide 3-kinase regulatory subunit 5 (PIK3R5)

<b>Product Code</b>	CSB-EP841225HU
<b>Relevance</b>	Regulatory subunit of the PI3K gamma complex. Required for recruitment of the catalytic subunit to the plasma mbrane via interaction with beta-gamma G protein dimers. Required for G protein-mediated activation of PIK3CG .
<b>Abbreviation</b>	Recombinant Human PIK3R5 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WYR1
<b>Alias</b>	PI3-kinase p101 subunitPhosphatidylinositol 4,5-bisphosphate 3-kinase regulatory subunit ;PtdIns-3-kinase regulatory subunitProtein FOAP-2;PtdIns-3-kinase p101p101-PI3K
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	MDSGYVEDSEESSSEWPWRRGSQERRGHRPGQKFIRIYKLFKSTSSQLVLRRL DSRSLEGSSDTALPLRRAGSLCSPLEPVSPPSRAQRSRSLPQPKLGTQLPS WLLAPASRPQRRRPFLSGDEDPKASTLRVVVFGSDRISGKVARAYSNLRRLE NNRPLLTRFFKLQFFYVPVKRSHGTSPGACPPPRSQTSPPTDSPRHASPGE LGTTPWEESTNDISHYLGMLDPWYERNVLGLMHLPPPEVLCQQSLKAEQAQALE GSPTQLPILADMLLYCRFAARPVLLQVYQTELTFTITGEKTTEIFIHSLELGHSAA TRAIKASGPGSKRLGIDGDREAVPLTLQIIYSKGAISGRSRWSNLEKVCTSVNL NKACRKQEELDSSMEALTLNLTEVVKRQNSKSKKGFNQISTSQIKVDKVIIGS NSCPFAVCLDQDERKILQSVVRCEVSPCYKPEKSDLSSPPQTPPDLPAAQAPD LCSSLCLPIMTFSGALP
<b>Research Area</b>	Signal Transduction
<b>Source</b>	E.coli
<b>Target Names</b>	PIK3R5
<b>Expression Region</b>	1-494aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	70.9kDa
<b>Protein Length</b>	Full Length of Isoform 2



## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.