



Recombinant Mouse Secretoglobin family 3A member 2 (Scgb3a2)

Product Code	CSB-MP846028MO
Abbreviation	Scgb3a2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q920H1
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	LLINRLPVV DKLPVPLDDI IPSFDPLKML LKTLGISVEH LVTGLKKCVD ELGPEASEAV KKLLVIIIICS YFPGRSLCYV NNLPSFVSVL FLPMICAYPR DSKKQTFAFI ERVFEQSKL
Source	Mammalian cell
Target Names	Scgb3a2
Protein Names	Recommended name: Secretoglobin family 3A member 2 Alternative name(s): Pneumo secretory protein 1 Short name= PnSP-1 Uteroglobin-related protein 1
Expression Region	22-139
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein of Isoform C
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.