





# Recombinant Escherichia coli O157:H7 Outer membrane protein C (ompC)


<b>Product Code</b>	CSB-EP848941EOD
<b>Relevance</b>	Forms pores that allow passive diffusion of small molecules across the outer membrane
<b>Abbreviation</b>	Recombinant E.coli O157:H7 ompC protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8XE41
<b>Alias</b>	Outer membrane protein 1B Porin OmpC
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O157:H7
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	AEVYNKDGKLDLYGKVDGLHYFSDDKSVDGDQTYMRLGFKGETQVTDQLT GYGQWEYQIQGNSAENENNSWTRVAFAGLKFQDVGSFDYGRNYGVVYDVTS WTDVLPFEGGDTYGSDNFMQQRGNGFATYRNTDFFGLVDGLNFAVQYQGKN GSVSGEGMTNNGREALRQNGDVGGSITYDYEGFGIGAAVSSSKRTDDQNS PLYIGNGDRAETYTGGLKYDANNIYLAAQYTQTYNATRVSGLGWANKAQNFE AVAQYQFDFGLRPSLAYLQSKGKNLGVINGRNYDDEDILKYVDVGATYYFNKN MSTYVDYKINLLDDNQFTRDAGINTDNIVALGLVYQF
<b>Research Area</b>	Microbiology
<b>Source</b>	E.coli
<b>Target Names</b>	ompC
<b>Protein Names</b>	Recommended name: Outer membrane protein C Alternative name(s): Outer membrane protein 1B Porin OmpC
<b>Expression Region</b>	22-367aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	58.4kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP848941EOD could indicate that this peptide derived from E.coli-expressed Escherichia coli O157:H7 ompC.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP848941EOD could indicate that this peptide derived from E.coli-expressed Escherichia coli O157:H7 ompC.

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.