



# Recombinant Mouse Cystathionine beta-synthase (Cbs)

<b>Product Code</b>	CSB-MP849665MO
<b>Abbreviation</b>	Cbs
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q91WT9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	PSGTSQCED GSAGGFQHLDMHSEKRQLEK GPSGDKDRVW IRPDTPSRCT WQLGRAMADS PHYHTVLTKS PKILPDILRK IGNTPMVRIN KISKNAGLKC ELLAKCEFFN AGGSVKDRIS LRMIEDAERA GNLKPGDTII EPTSGNTGIG LALAAVKG Y RCIIVMPEKM SMEKVDV LRA LGAEIVR TPT NARFDSPESH VGVAWRLKNE IPNSHILDQY RNASNPLAHY DDTAE EILQQ CDGKLDMLVA SAGTGGTITG IARKLKEKCP GCKIIGVDPE GSILAEPEEL NQTEQTAYEV EGIGYDFIPT VLDRAVVDKW FKSND EDSFA FARMLIAQEG LLCGGSSGSA MAVAVKAARE LQEGQR CVVI LPDSVR NYMS KFLSDKWMLQ KGFMKEELSV KRPWWWRLRV QELSLSAPLT VLPTVTCEDT IAILREK GFD QAPVVNESGA ILGMVTLGNM LSSLLAGKVR PSDEVCKVLY KQFKPIHLTD TLGTLSHILE MDHFALVVHE QIQSRDQAWS GVVGGPTDCS NGMSSKQQMV FGVVTAIDLL NFVAAREQTQ T
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Cbs
<b>Protein Names</b>	Recommended name: Cystathionine beta-synthase EC= 4.2.1.22 Alternative name(s): Beta-thionase Serine sulfhydrase
<b>Expression Region</b>	2-561
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein acts as a homotetramer to catalyze the conversion of homocysteine to cystathionine, the first step in the transsulfuration pathway. The encoded protein is allosterically activated by adenosyl-methionine and uses pyridoxal phosphate as a cofactor. Defects in this gene can cause cystathionine beta-synthase deficiency (CBS), which can lead to homocystinuria.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.