



# Recombinant Mouse Alanine aminotransferase 1 (Gpt)

<b>Product Code</b>	CSB-EP851235MO
<b>Abbreviation</b>	Gpt
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8QZR5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ASQRNDRIQ ASRNGCLKGKV LTLDTMNPCV RRVEYAVRGP IVQRALELEQ ELRQGVKKPF TEVIRANIGD AQAMGQRPIT FFRQVLALCV YPNLLSSPDF PEDAKRRAER ILQACGGHSL GAYSISSGIQ PIREVAQYI ERRDGGIPAD PNNIFLSTGA SDAIVTMLKL LVAGEGRART GVLIPQYP LYSAAEAELD AVQVDYYLDE ERAWALDIAE LRRALCQARD RCCPRVLCVI NPGNPTGQVQ TRECIEAVIR FAFEEGLFLM ADEVYQDNVY AEGSQFHSFK KVLTEMGPPY ATQQELASFH SVSKGYMGEC GFRGGYVEVV NMDAEVQKQM AKLMSVRLCP PVPGQALMGM VVSPPTPSEP SFKQFQAERQ EVLAELAACA KLTEQVFNEA PGIRCNPVQG AMYSFPQIQL PLKAVQRAQD LGLAPDMFFC LCLLEETGIC VVPGSGFGQQ EGTYHFRMTI LPPMEKLRVL LEKLRHFHAK FTHEYS
<b>Source</b>	E.coli
<b>Target Names</b>	Gpt
<b>Protein Names</b>	Recommended name: Alanine aminotransferase 1 Short name= ALT1 EC= 2.6.1.2 Alternative name(s): Glutamate pyruvate transaminase 1 Short name= GPT 1 Glutamic--alanine transaminase 1 Glutamic--pyruvic transaminase 1
<b>Expression Region</b>	2-496
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.