



Recombinant Mouse Alanine aminotransferase 1 (Gpt)

Product Code	CSB-MP851235MO
Abbreviation	Gpt
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8QZR5
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	ASQRNDRIQ ASRNGLGKGV LTLDTMNPCV RRVEYAVRGP IVQRALELEQ ELRQGVKKPF TEVIRANIGD AQAMGQRPIT FFRQVLALCV YPNLLSSPDF PEDAKRRAER ILQACGGHSL GAYSISSGIQ PIREVAQYI ERRDGGIPAD PNNIFLSTGA SDAIVTMLKL LVAGEGRART GVLIPQYP LYSAAEAELD AVQVDYYLDE ERAWALDIAE LRRALCQARD RCCPRVLCVI NPGNPTGQVQ TRECIEAVIR FAFEEGLFLM ADEVYQDNVY AEGSQFHSFK KVLTEMGPPY ATQQELASFH SVSKGYMGEC GFRGGYVEVV NMDAEVQKQM AKLMSVRLCP PVPGQALMGM VVSPPTPSEP SFKQFQAERQ EVLAELAACA KLTEQVFNEA PGIRCNPVQG AMYSFPQIQL PLKAVQRAQD LGLAPDMFFC LCLLEETGIC VVPGSGFGQQ EGTYHFRMTI LPPMEKLRVL LEKLRHFHAK FTHEYS
Source	Mammalian cell
Target Names	Gpt
Protein Names	Recommended name: Alanine aminotransferase 1 Short name= ALT1 EC= 2.6.1.2 Alternative name(s): Glutamate pyruvate transaminase 1 Short name= GPT 1 Glutamic--alanine transaminase 1 Glutamic--pyruvic transaminase 1
Expression Region	2-496
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.