



Recombinant Human SARS coronavirus Replicase polyprotein 1ab (rep), partial, Yeast

Product Code	CSB-YP314477HQE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C6X7
Product Type	Recombinant Protein
Immunogen Species	Human SARS coronavirus (SARS-CoV) (Severe acute respiratory syndrome coronavirus)
Purity	≥85% (SDS-PAGE)
Source	Yeast
Target Names	rep
Protein Names	Recommended name: Replicase polyprotein 1ab Short name= pp1ab Alternative name(s): ORF1ab polyprotein Cleaved into the following 15 chains: 1. Non-structural protein 1 Short name= 2. nsp1 Alternative name(s): Leader protein Non-str
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged Tag-Free
Protein Length	partial, Nsp13 function region
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.