



Recombinant Alkaliphilus metalliredigens Glutaminase (glsA)

Product Code	CSB-YP009528AUP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6TU96
Product Type	Recombinant Protein
Immunogen Species	Alkaliphilus metalliredigens (strain QYMF)
Purity	>85% (SDS-PAGE)
Sequence	MNEQLQKILD TNRHIIQGGQ LPTYIPELSK ANKEALGIYV ADLDGNEYGV GDYEYPFTIQ SISKVVTLLL ALSDRGEKYV FDKVGMEPTG DPFNSMMKLE VVRPSKPFNP MINAGIAVT SMIKGDSQKE RLERILGFFR QLTENPNLQV NQSVYRSEKI TGDRNRSMAY FMKDVGIMQN DIESDLDL YF NQCSIEVTCK DIAKIGRFLA NGGVLFETGQ QLIDEKYIKM AEA FMVTCGM YNASGEFAIK VGIPAKSGVG GGIMASVPPK MGIGVVGPSL DEKGN SIAGV RVLQKMAEDY GLSIF
Source	Yeast
Target Names	glsA
Protein Names	Recommended name: Glutaminase EC= 3.5.1.2
Expression Region	1-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.