



Recombinant Human Store-operated calcium entry-associated regulatory factor (TMEM66)

Product Code	CSB-EP853392HU-B
Abbreviation	TMEM66
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96BY9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	SDGQYS PPPYSEYPPF SHRYQRFTNS AGPPPPGFKS EFTGPQNTGH GATSGFGSAF TGQQGYENSG PGFWTGLGTG GILGYLFGSN RAATPFSDSW YPSYPPSYG GTWNRAYSPL HGGSGSYSVC SNSDTKTRTA SGYGGTRRR
Source	E.coli
Target Names	SARAF
Protein Names	Recommended name: Store-operated calcium entry-associated regulatory factor Short name= SARAF Short name= SOCE-associated regulatory factor Alternative name(s): HBV X-transactivated gene 3 protein HBV XAg-transactivated protein 3
Expression Region	195-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.