




Recombinant Human Protein boule-like (BOLL)

Product Code	CSB-EP854082HU
Relevance	Probable RNA-binding protein, which may be required during spermatogenesis. May act by binding to the 3'-UTR of mRNAs and regulating their translation
Abbreviation	Recombinant Human BOLL protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N9W6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	<p>MQTDSLSPSPNPVSPVPLNNTSAPRYGTVIPNRIFVGGIDFKTNESDLRKFFS QYGSVKEVKIVNDRAGVSKGYGFVTFETQEDAQKILQEAELNYKDKKLNIGP AIRKQQVGIPRSSIMPAAGTMYLTTSTGYPTYHNGVAYFHTPEVTSVPPPWP SRVCSPPVMVAQPIYQQPAYHYQATTQYLPQQWQWSVPQPSASSAPFLYL QPSEVIYQPVEIAQDGGCVPPPLSLMETSVPPEPYS DHGVQATYHQVYAPSAIT MPAPVMQPEPIKTVWSIH</p>
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	BOLL
Protein Names	Recommended name: Protein boule-like
Expression Region	1-283aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	58.3kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.