




Recombinant Human Coiled-coil domain-containing protein 3 (CCDC3)

Product Code	CSB-EP863099HU
Abbreviation	Recombinant Human CCDC3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BQI4
Alias	Fat/vessel-derived secretory protein Short name: Favine
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	CQLPSEWRPLSEGCRAELAETIVYARVLALHPEAPGLYNHLPWQYHAGQGGL FYSAEVEMLCDQAWGSMLEVPAGSRLNLTGLGYFSCHSHTVVQDYSYFFFLR MDENYNLLPHGVNFQDAIFPDTQENRRMFSSLFQFSNCSQGQQLATFSSDW EIQEDSRLMCSSVQKALFEEEDHVKKLQQKVATLEKRNRLRERVKKVKRSLR QARKKGRHLELANQKLSEKLAAGALPHINARGPVRPPYLRG
Research Area	Others
Source	E.coli
Target Names	CCDC3
Protein Names	Recommended name: Coiled-coil domain-containing protein 3 Alternative name(s): Fat/vessel-derived secretory protein Short name= Favine
Expression Region	22-270aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	32.1kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.