




Recombinant Mouse Interleukin-23 subunit alpha (IL23a)

Product Code	CSB-EP863641MO
Relevance	Associates with IL12B to form the IL-23 interleukin, a heterodimeric cytokine which functions in innate and adaptive immunity. IL-23 may constitute with IL-17 an acute response to infection in peripheral tissues. IL-23 binds to a heterodimeric receptor complex composed of IL12RB1 and IL23R, activates the Jak-Stat signaling cascade, stimulates memory rather than naive T-cells and promotes production of proinflammatory cytokines. IL-23 induces autoimmune inflammation and thus may be responsible for autoimmune inflammatory diseases and may be important for tumorigenesis.
Abbreviation	IL23a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9EQ14
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	VPRSSSPDWAQCQQLSRNLCMLAWNAHAPAGHMNLLREEEDEETKNNVPRI QCEDGCDPQGLKDNSQFCLQRIRQGLAFYKHLLEDSDIFKGEALLPDSPMEQL HTSLLGLSQLLQPEDHPRETQQMPSLSSSQWQRPLLRSKILRSLQAFLAIAA RVFAHGAATLTEPLVPTA
Research Area	Immunology
Source	E.coli
Target Names	IL23a
Protein Names	Recommended name: Interleukin-23 subunit alpha Short name= IL-23 subunit alpha Short name= IL-23-AA Alternative name(s): Interleukin-23 subunit p19 Short name= IL-23p19
Expression Region	22-196
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	23.7kDa
Protein Length	Full length of mature protein



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.