



Recombinant Human L-2-hydroxyglutarate dehydrogenase, mitochondrial (L2HGDH)

Product Code	CSB-EP864008HU
Abbreviation	Recombinant Human L2HGDH protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H9P8
Alias	Duranin
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	VIVGGGIVGLASARALILRHPSLSIGVLEKEKDLAVHQTGHNSGVIHSGIYYKPE SLKAKLCVQGAALLYEYCQQKGISYKQCGKLIVAVEQEEIPRLQALYEKGLQNG VPGLRLIQQEDIKKKEPYCRGLMAIDCPHTGIVDYRQVALSFAQDFQEAGGSVL TNFEVKGIEMAKESPSRSIDGMQYPIVIKNTKGEEIRCQYVVT CAGLYSDRISEL SGCTPDPRIVPFRGDYLLKPEKCYLVKGNIYPVPSRFPFLGVHFTPRMDGSI WLGPNVAVLAFKREGYRPFDFSATDVMIIINSGLIKLASQNFSGVTEMYKACF LGATVKYLQKFIPEITISDILRGPAGVRAQALDRDGNLVEDFVFDAGVGDIGNRI LHVRNAPSPAATSSIAISGMIADDEVQQR FEL
Research Area	Cell Biology
Source	E.coli
Target Names	L2HGDH
Protein Names	Recommended name: L-2-hydroxyglutarate dehydrogenase, mitochondrial EC=1.1.99.2 Alternative name(s): Duranin
Expression Region	52-463aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	61.3kDa
Protein Length	Full Length of Mature Protein
Image	

