



Recombinant Human E3 ubiquitin-protein ligase TRAIP (TRAIP)

Product Code	CSB-EP866296HU
Relevance	E3 ubiquitin ligase acting as a negative regulator of innate immune signaling. Inhibits activation of NF-kappa-B mediated by TNF. Negatively regulates TLR3/4- and RIG-I-mediated IRF3 activation and subsequent IFNB1 production and cellular antiviral response by promoting 'Lys-48'-linked polyubiquitination of TNK1 leading to its proteasomal degradation. Involved in response to genotoxic lesions during genome replication. Promotes H2AX and RPA2 phosphorylation after replication-associated DNA damage and assists fork progression at UV-induced replication-blocking lesions during S phase. Has also been proposed to play a role in promoting translesion synthesis by mediating the assembly of 'Lys-63'-linked poly-ubiquitin chains on the Y-family polymerase POLN in order to facilitate bypass of DNA lesions and preserve genomic integrity. The function in translesion synthesis is controversial
Abbreviation	Recombinant Human TRAIP protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BWF2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MPIRALCTICSDFFDHSRDVAAIHCGHTFHLQCLIQWFETAPSRTCPQCRIQVG KRTIINKLFFDLAQEEENVLDAEFLKNELDNVRAQLSQKDKEKRDSQVIIDTLRD TLEERNATVVSLQQALGKAEMLCSTLKKQMKYLEQQQDETQQAQEEARRLRS KMKTMEQIELLLQSQRPEVEEMIRDMGVGQSAVEQLAVYCVSLKKEYENLKEA RKASGEVADKLRKDLFSSRSKLQTVYSELDQAKLELKSAQKDLQSADKEIMSL KKKLTMLQETLNLPPVASETVDRVLVLESPAPVEVNLKLRPSFRDDIDLNATFD VDTPPARPSSSQHGYYEKLCKLEKSHSPIQDVPKKICKGPRKESQLSLGGQSCA GEPDEELVGAFPIFVRNAILGQKQPKRPRSESSCSKDVVRTGFDGLGGRTKFI QPTDTVMIRPLPVKPKTKVKQVRVKTVPVSLFQAKLDTFLWS
Research Area	Others
Source	E.coli
Target Names	TRAIP
Protein Names	RING finger protein 206 RING-type E3 ubiquitin transferase TRAIPCurated TRAF-interacting protein RNF206, TRIP
Expression Region	1-469aa



Notes

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Tag Info

N-terminal 10xHis-tagged


Mol. Weight

58.8 kDa


Protein Length

Full Length


Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP866296HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) TRAIP.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP866296HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) TRAIP.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.