



Recombinant Dog Bone morphogenetic protein 7 (BMP7)

Product Code	CSB-MP002744DO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34819
Product Type	Recombinant Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Purity	≥85% (SDS-PAGE)
Sequence	GKHNSAPMFM LDLYNAMAVE EGGGPAGQGF SYPYKAVFST QGPPLASLQD SHFLTADAMV MSFVNLEHD KEFFHPRYHH REFRFDLSKI PEGEAVTAAE FRIYKDYIRE RFDNETFRIS VYQVLQEH LG RESDLFLLDS RTLWASEEGW LVFDITATSN HWVVNPRHNL GLQLCVETLD GQSINPK
Source	Mammalian cell
Target Names	BMP7
Protein Names	Recommended name: Bone morphogenetic protein 7 Short name= BMP-7 Alternative name(s): Osteogenic protein 1 Short name= OP-1
Expression Region	1-187
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	The bone morphogenetic proteins (BMPs) are a family of secreted signaling molecules that can induce ectopic bone growth. Many BMPs are part of the transforming growth factor-beta (TGFB) superfamily. BMPs were originally identified by an ability of demineralized bone extract to induce endochondral osteogenesis in vivo in an extraskeletal site. Based on its expression early in embryogenesis, the BMP encoded by this gene has a proposed role in early development and possible bone inductive activity.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.