



# Recombinant Human Gamma-1-syntrophin (SNTG1)

<b>Product Code</b>	CSB-MP873648HU
<b>Abbreviation</b>	SNTG1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NSN8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	MDFRTACEET KTGICLLQDG NQEPFKVRLH LAKDILMIQE QDVICVSGEP FYSGERTVTI RRQTVGGFGL SIKGGAEHNI PVVVSKE QRAELSGLLF IGDAILQING INVRKCRHEE VVQVLRNAGE EVTLTVSFLK RAPAFKLPL NEDCACAPSD QSSGTSSPLC DSGHLHLYHP NNTDTLSCSS WPTSPGLRWE KRWCDLRLIP LLHSRFSQYV PGTDLSRQNA FQVIAVDGVC TGIIQCLSAE DCVDWLQAIA TNISNLTKHN IKKINRNFV NQQIVYMGWC EAREQDPLQD RVYSPTFLAL RGSCLYKFLA PPVTTWDWTR AEKTFVSYEI MCKILKDSL LDRRKQCFTV QSEGEDLYF SVELESLAQ WERAFQTATF LEVERIQCKT YACVLESHLM GLTIDFSTGF ICFDAATKAV LWRYKFSQLK GSSDDGKSKI KFLFQNPDTK QIEAKELEFS NLFVHLCHH SFFAAKVAQL DPLFLGNQAT ASTAASSATT SKAKYTT
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SNTG1
<b>Protein Names</b>	Recommended name: Gamma-1-syntrophin Short name= G1SYN Alternative name(s): Syntrophin-4 Short name= SYN4
<b>Expression Region</b>	1-517
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.