



Recombinant Staphylococcus epidermidis Endoribonuclease MazF (mazF)

Product Code	CSB-EP880696FLL
Abbreviation	Recombinant Staphylococcus epidermidis mazF protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9F7V5
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Staphylococcus epidermidis
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MIRRGDVYLADLSPVQQGSEQGGVRPVVIIQNNTGNKYSPVIVAAITDGINKAKI PTHVEIEKKKYKLDKDSVILLEQIRTLDDKRLKEKLTFLSESKMIEVDNALDISLG LNNFDHHKS
Research Area	Others
Source	E.coli
Target Names	mazF
Protein Names	Recommended name: mRNA interferase MazF EC= 3.1.-.-Alternative name(s): Endoribonuclease MazF Toxin MazF
Expression Region	1-120aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	19.0 kDa
Protein Length	Full Length
Image	

