




# Recombinant Human 2-Hydroxyacid oxidase 2 (HAO2)

<b>Product Code</b>	CSB-EP882145HU
<b>Relevance</b>	Catalyzes the oxidation of L-alpha-hydroxy acids as well as, more slowly, that of L-alpha-amino acids.
<b>Abbreviation</b>	Recombinant Human HAO2 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NYQ3
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	SLVCLTDFQAHAREQLSKSTRDFIEGGADDSITRDDNIAAFKRIRLRPRYL RDV SEVDTRTTIQGEEISAPICIAPTGFHCLVWPDGEMSTARAAQAAGICYITSTFAS CSLEDIVIAAPEGLRWFQLYVHPDLQLNKQLIQRVESLGFKALVITLDTVPVCGNR RHDIRNQLRRNLTLTDLQSPKKGNAIPYFQMTPISTSLCWNDLSWFQSITRLPII LKGILTKEAELAVKHNVQGIIVSNHGGQRQLDEVLASIDALTEVVAAVKKGKIEVY LDGGVRTGNDVLKALALGAKCIFLGRPILWGLACKGEHGVKEVLNILTNEFHTS MALTGCRSVAEINRNLVQFSRL
<b>Research Area</b>	Tags & Cell Markers
<b>Source</b>	E.coli
<b>Target Names</b>	HAO2
<b>Protein Names</b>	Recommended name: Hydroxyacid oxidase 2 Short name= HAOX2 EC= 1.1.3.15 Alternative name(s): (S)-2-hydroxy-acid oxidase, peroxisomal Cell growth-inhibiting gene 16 protein Long chain alpha-hydroxy acid oxidase Long-chain L
<b>Expression Region</b>	2-351aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	54.7kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.