



# Recombinant *Alkaliphilus oremlandii* Uracil-DNA glycosylase (ung)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP025641AUQ   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | A8MM35  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Alkaliphilus oremlandii</i> (strain OhILAs) ( <i>Clostridium oremlandii</i> (strain OhILAs))   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MGVIFENDWQ ELLKDEFEKE YYKKVRGFLV EEYRTKTIYP DKNHIFSALH<br>YTAYKDVRV V ILGQDPYHGP NQAHGLSFSV QPGVRIPPSL VNIYKELET<br>D LGCYIPNHGY LKKWADQGV L LLNTSLTVRA GEANSHQKIG WEIFTDHIIQ<br>LLNDRKEPIV FILWGNN AIS KERFITNPQH YIIKSVHPSP LSASRGFFGS<br>KPF SKTNSFL ESIGEKPIDW QIENR  |
| <b>Source</b>            | <i>E. coli</i>  |
| <b>Target Names</b>      | ung   |
| <b>Protein Names</b>     | Recommended name: Uracil-DNA glycosylase Short name= UDG EC= 3.2.2.27   |
| <b>Expression Region</b> | 1-225   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |