




Recombinant Human SAP30-binding protein (SAP30BP)

Product Code	CSB-EP883403HU
Relevance	Induces cell death. May act as a transcriptional corepressor of a gene related to cell survival. May be involved in the regulation of beta-2-microglobulin genes.
Abbreviation	SAP30BP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UHR5
Alias	Transcriptional regulator protein HCNGP
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MAGKKNVLSSLAVYAEDSEPESEDGEAGIEAVGSAAEEKGGLVSDAYGEDDFS RLGGDEDGYEEEEEDENSRQSEDDDDSETEKPEADDPKDNTAEKRDPQELVA SFSERVRNMPDEIKIPPEPPGRCSNHLQDKIQKLYERKIKEGMMDMNYIIQRKK EFRNPSIYEKLIQFCAIDELGTNYPKDMFDPHGWSEDSYYEALAKAQKIEMDKL EKAKKERTKIEFVTGTTKGGTTTNTATSTTTTTASTAVADAQKRKSKWDSAIPVTTI AQPTILTTTATLPAVVTVTTTSASGSKTTVISAVGTIVKKAKQ
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	SAP30BP
Protein Names	Recommended name: SAP30-binding protein Alternative name(s): Transcriptional regulator protein HCNGP
Expression Region	1-308aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	60.9kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.