




Recombinant Mouse Thioredoxin domain-containing protein 12 (Txndc12)

Product Code	CSB-EP887432MO
Relevance	Possesses significant protein thiol-disulfide oxidase activity.
Abbreviation	Recombinant Mouse Txndc12 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9CQU0
Alias	Endoplasmic reticulum resident protein 19 Short name: ER protein 19 Short name: ERp19 Thioredoxin-like protein p19
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	RTGLGKGFGDHIHWRTLEDGKKKAAASGLPLMVIIHKSWCGACKALKPKFAES TEISELSHNFVMVNLEDEEPRDEDFSPDGGYIPRILFLDPSGKVRPEIINESGN PSYKYFYVSAEQVVQGMKEAQRRLTGDAFREKHFQDEL
Research Area	Others
Source	E.coli
Target Names	Txndc12
Expression Region	25-170aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	20.5kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.