



Recombinant Pig Lactadherin (MFGE8)

Product Code	CSB-MP013752PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P79385
Product Type	Recombinant Protein
Immunogen Species	Sus scrofa (Pig)
Purity	>85% (SDS-PAGE)
Sequence	FSGDFCDSSL CLNGGTCLLD QDPQKPFHCL CPEGFTGLIC NETEKGPCFP NPCHNDAECE VIDDAHRGDV FTEYICKCPH GYTGIHCEII CNAPLGMETG AIADFAQISAS SMHLGFMGLQ RWAPELARLH RAGIVNAWTA SNYDRNPWIQ VNULLRRMRVT GVVVTQGASRA GSAEYMKTFK VAYSTDGRKF QFIQGAEESG DKIFMGNLDN SGLKVNLFEEV PLEVQYVRLV PIICHRGCTL RFELLGCELS GCAEPLGLKD NTIPNKQITA SSFYRTWGLS AFSWYPFYAR LDNQGKFNAW TAQNSASEW LQIDLGSQRR VTGIITQGAR DFGHIQYVAA YKVAYSDDGV SWTEYRDQGA LEGKIFPGNL DNNSHKKNMF ETPFLTRFVR ILPVAWHNRI TLRVELLGC
Source	Mammalian cell
Target Names	MFGE8
Protein Names	Recommended name: Lactadherin Alternative name(s): MFGM Milk fat globule-EGF factor 8 Short name= MFG-E8 SED1 Sperm surface protein SP47 Short name= PP47
Expression Region	1-409
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.