




Recombinant Human Cysteine-rich hydrophobic domain 2 protein (CHIC2)

Product Code	CSB-EP892453HU
Abbreviation	Recombinant Human CHIC2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UKJ5
Alias	BrX-like translocated in leukemia
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MADFDEIYEEEEDEERALEEQLLKYSPPVVRGSGHVTVFGLSNKFESEFPS SLTGKVAPEEFKASINRVNSCLKKNLPVNVRLWLLCGCLCCCCTLGCSMWPVIC LSKRTRRSIEKLLEWENNRLYHKLCLHWRLSKRKCETNNMMEYVILIEFLPKTPI FRPD
Research Area	others
Source	E.coli
Target Names	CHIC2
Protein Names	Recommended name: Cysteine-rich hydrophobic domain 2 protein Alternative name(s): BrX-like translocated in leukemia
Expression Region	1-165aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	46.3kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.