



Recombinant Human CD9 antigen (CD9), partial

Product Code	CSB-BP004969HU1
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P21926
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SHKDEVIKEVQEFYKDTYNKCLKTKDEPQRETLKAIHYALNCCGLAGGVEQFISD ICPKKD VLETFTVKSCPDAIKEVFDNKFHI
Source	Baculovirus
Target Names	CD9
Protein Names	Recommended name: CD9 antigen Alternative name(s): 5H9 antigen Cell growth-inhibiting gene 2 protein Leukocyte antigen MIC3 Motility-related protein Short name= MRP-1 Tetraspanin-29 Short name= Tspan-29 p24 CD_antigen=
Expression Region	112-195aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.