




Recombinant Mouse Acid ceramidase (Asah1), partial

Product Code	CSB-EP895296MO
Relevance	Hydrolyzes the sphingolipid ceramide into sphingosine and free fatty acid.
Abbreviation	Recombinant Mouse Asah1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9WV54
Alias	Acylsphingosine deacylase N-acylsphingosine amidohydrolase
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	QAQDVPPWTEDCRKSTYPPSGPTYRGPVPWHTINLDLPPYKRWHELLAQKAP ALRILVNSITSLVNTFVPSGKLMKMVDQKLPGMIGSLPDPFGEEEMRGIADVTGI PLGEIISFNIFYELFTM
Research Area	Signal Transduction
Source	E.coli
Target Names	Asah1
Protein Names	Recommended name: Acid ceramidase Short name= AC Short name= ACDase Short name= Acid CDase EC= 3.5.1.23 Alternative name(s): Acylsphingosine deacylase N-acylsphingosine amidohydrolase Cleaved into the following 2 ch
Expression Region	19-141aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10XHis-GST-tagged and C-terminal Myc-tagged
Mol. Weight	43.8kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.