




# Recombinant Mouse CCN family member 5 (Ccn5)

|                          |                                                                                                                                                                                                                                                                                 |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP896674MO                                                                                                                                                                                                                                                                  |
| <b>Relevance</b>         | May play an important role in modulating bone turnover. Promotes the adhesion of osteoblast cells and inhibits the binding of fibrinogen to integrin receptors. In addition, inhibits osteocalcin production .                                                                  |
| <b>Abbreviation</b>      | Recombinant Mouse Ccn5 protein                                                                                                                                                                                                                                                  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| <b>Uniprot No.</b>       | Q9Z0G4                                                                                                                                                                                                                                                                          |
| <b>Alias</b>             | CCN family member 5Connective tissue growth factor-like protein ;CTGF-L                                                                                                                                                                                                         |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                             |
| <b>Immunogen Species</b> | Mus musculus (Mouse)                                                                                                                                                                                                                                                            |
| <b>Purity</b>            | ≥ 90% as determined by SDS-PAGE.                                                                                                                                                                                                                                                |
| <b>Sequence</b>          | QLCPAPCACPWTPPPQCPPGVPLVLDGCGCCRVLCARRLGESCDHLHVCDPSQ<br>GLVCQPGAGPSGRGAVCLFEEDDGSCEVNGRRYLDGETFKPNCRVLCRCDD<br>GGFTCLPLCEDVRLPSWDCPRPRRIQVPGRCCEWVCDQAVMQPAIQPSS<br>AQGHQLSALVTPASADGPCPNWSTAWGPCSTTCGLGIATRVSQNQRFCQLEI<br>QRRCLLSRPCLASRSHGSWNSAF                            |
| <b>Research Area</b>     | Others                                                                                                                                                                                                                                                                          |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                          |
| <b>Target Names</b>      | Ccn5                                                                                                                                                                                                                                                                            |
| <b>Protein Names</b>     | Recommended name: WNT1-inducible-signaling pathway protein 2 Short name= WISP-2Alternative name(s): CCN family member 5 Connective tissue growth factor-like protein Short name= CTGF-L                                                                                         |
| <b>Expression Region</b> | 24-251aa                                                                                                                                                                                                                                                                        |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                             |
| <b>Tag Info</b>          | N-terminal 6xHis-SUMO-tagged                                                                                                                                                                                                                                                    |
| <b>Mol. Weight</b>       | 40.5kDa                                                                                                                                                                                                                                                                         |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                   |
| <b>Image</b>             |                                                                                                                                                                                                                                                                                 |



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.