



Recombinant Human Receptor-interacting serine/threonine-protein kinase 3 (RIPK3)

Product Code	CSB-BP897497HU
Relevance	Essential for necroptosis, a programmed cell death process in response to death-inducing TNF-alpha family members. Upon induction of necrosis, RIPK3 interacts with, and phosphorylates RIPK1 and MLKL to form a necrosis-inducing complex. RIPK3 binds to and enhances the activity of three metabolic enzymes: GLUL, GLUD1, and PYGL. These metabolic enzymes may eventually stimulate the tricarboxylic acid cycle and oxidative phosphorylation, which could result in enhanced ROS production.
Abbreviation	Recombinant Human RIPK3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y572
Alias	RIP-like protein kinase 3 Receptor-interacting protein 3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MSCVKLWPSGAPAPLVSIIELENQELVGKGGFGTVFRAQHRKWGYDVAVKIV NSKAISREVKAMASLDNEFVLRLEGVIEKVNWDQDPKALVTKFMENGSLSG LQSQCPRPWPLLCRLLKEVVLGMFYLDQNPVLLHRDLKPSNVLLDPELHVKL ADFGGLSTFQGGSSQSGTGSSEPGGTGLGYLAPELFVNVNRKASTASDVYSFGIL MWAVLAGREVELPTEPSLVYEAVCNQRNRPDLAELPQAGPETPGLEGLKELM QLCWSSEPKDRPSFQECLPKTDEVFQMVENNMNAAVSTVKDFLSQLRSSNR RFSIPESGQGGTEMDGFRRTIENQHRSRNDVMVSEWLNKLNLEPPSSVPKKC PSLTKRSRAQEEQVPQAWTAGTSSDSMAQPPQTPETSTFRNQMPSPSTSTGT PSPGPRGNQGAERQGMNWSCRTPEPNPVTGRPLVNIYNCSGVQVGDNNYLT MQQTALPTWGLAPSGKGRGLQHPPVGSQEGPKDPEAWSRPQGWYNHS GK
Research Area	Cell Biology
Source	Baculovirus
Target Names	RIPK3
Protein Names	Recommended name: Receptor-interacting serine/threonine-protein kinase 3 EC= 2.7.11.1Alternative name(s): RIP-like protein kinase 3 Receptor-interacting protein 3 Short name= RIP-3
Expression Region	1-518aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at



4°C for up to one week.

Tag Info


N-terminal 10xHis-tagged and C-terminal Myc-tagged

Mol. Weight

60.9kDa

Protein Length

Full Length

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.