




Recombinant Streptococcus pneumoniae serotype 4 Pneumolysin (ply)


Product Code	CSB-BP314690FMW
Relevance	Sulfhydryl-activated toxin that causes cytolysis by forming pores in cholesterol containing host membranes. After binding to target membranes, the protein undergoes a major conformation change, leading to its insertion in the host membrane and formation of an oligomeric pore complex. Cholesterol may be required for binding to host membranes, membrane insertion and pore formation. Can be reversibly inactivated by oxidation
Abbreviation	Recombinant Streptococcus pneumoniae serotype 4 ply protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C2J9
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pneumoniae serotype 4 (strain ATCC BAA-334 / TIGR4)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	ANKAVNDFILAMNYDKKKLLTHQGESIENRFIKEGNQLPDEFVVIERKKRSLST NTSDISVTATNDSRLYPGALLVDETLLENNPTLLAVDRAPMTYSIDLPLGLASS DSFLQVEDPSNSSVRGAVNDLLAKWHQDYGQVNNVPARMQYEKITAHSMEQ LKVKFGSDFEKTGNSLDIDFNSVHSGEKQIQIVNFKQIYYTVSVDVKNPGDVF QDVTVEDLKQRGISAERPLVYISSVAYGRQVYLKLETTSKSDEVEAAFEALIK GVKVAPQTEWKQILDNTEVKAVILGGDPSSGARVVTGKVDMMVEDLIQEGSRFT ADHPGLPISYTTSSFLRDNVVATFQNSTDYVETKVTAYRNGDLLLDHSGAYVAQ YYITWNELSYDHQGKEVLTPKAWDRNGQDLTAHFTTSIPLKGNVRNLSVKIRE CTGLAWEWVRTVYEKTDLPLVRKRTISIWGTTLYPQVEDKVEND
Research Area	Others
Source	Baculovirus
Target Names	ply
Protein Names	Thiol-activated cytolysin
Expression Region	2-471aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	56.8 kDa
Protein Length	Full Length of Mature Protein




Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of Baculovirus host and target protein, the LC-MS/MS Analysis result of CSB-BP314690FMW could indicate that this peptide derived from Baculovirus-expressed Streptococcus pneumoniae serotype 4 (strain ATCC BAA-334 / TIGR4) ply.



Based on the SEQUEST from database of Baculovirus host and target protein, the LC-MS/MS Analysis result of CSB-BP314690FMW could indicate that this peptide derived from Baculovirus-expressed Streptococcus pneumoniae serotype 4 (strain ATCC BAA-334 / TIGR4) ply.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.