





Recombinant Human Cancer/testis antigen 1 (CTAG1A)

Product Code	CSB-EP006125HU
Abbreviation	Recombinant Human CTAG1A protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P78358
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MQAEGRGTTGGSTGDADGPGGPGIPDGPGGNAGGPGEAGATGGRGPRGAG AARASGPGGGAPRGPHGGAASGLNGCCRCGARGPESRLLEFYLAMPFATPM EAELARRSLAQDAPPLPVPGVLLKEFTVSGNILTIRLTAADHRQLQLSISSCLQQ LSLLMWITQCFLPVFLAQPPSGQRR
Research Area	others
Source	E.coli
Target Names	CTAG1A
Protein Names	Autoimmunogenic cancer/testis antigen NY-ESO-1Cancer/testis antigen 6.1
Expression Region	1-180aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	21.6kDa
Protein Length	Full Length
Image	



The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.