



Recombinant Human Low-density lipoprotein receptor-related protein 2 (LRP2), partial

Product Code	CSB-BP013096HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P98164
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	NCTASQFKCASGDKCIGVTNRCDGVFDCSDNSDEAGCPTRPPGMCHSDEFQ CQEDGICIPNFWECDGHPDCLYGSDEHNACVPKTCPPSSYFHCDNGNCIHRAW LCDRDNDCGDMSDEKDCPTQPFRCPQSWQWQCLGHNICVNLSVVCDFGIFDCP NGTDESPLCNGNSCDFNGGCTHECVQEPFGAKCLCPLGFLLANDSKTCE
Research Area	Immunology
Source	Baculovirus
Target Names	LRP2
Protein Names	Recommended name: Low-density lipoprotein receptor-related protein 2 Short name= LRP-2 Alternative name(s): Glycoprotein 330 Short name= gp330 Megalin
Expression Region	1186-1389AA
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.