



Recombinant Human Granzyme B (GZMB)

Product Code	CSB-MP010082HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P10144
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	IIGGHEAKPHSRPYMAYLMIWDQKSLKRCGGFLIRDDFVLTAAHCWGSSINVTL GAHNIKEQEPTQQFIPVKRPIHPAYNPKNFSNDIMLLQLERKAKRTRAVQPLR LPSNKAQVKPGQTCSVAGWGQTAPLGKHSHTLQEVKMTVQEDRKCESDLRH YYDSTIELCVGDPEIKKTSFKGDSGGPLVCNKVAQGIVSYGRNNGMPRACTK VSSFVHWIKKTMKRY
Research Area	Immunology
Source	Mammalian cell
Target Names	GZMB
Protein Names	Recommended name: Granzyme B EC= 3.4.21.79 Alternative name(s): C11 CTLA-1 Cathepsin G-like 1 Short name= CTSSL1 Cytotoxic T-lymphocyte proteinase 2 Short name= Lymphocyte protease Fragmentin-2 Granzyme-2
Expression Region	21-247aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Cytolytic T lymphocytes (CTL) and natural killer (NK) cells share the remarkable ability to recognize, bind, and lyse specific target cells. They are thought to protect their host by lysing cells bearing on their surface nonself antigens, usually peptides or proteins resulting from infection by intracellular pathogens. This protein is crucial for the rapid induction of target cell apoptosis by CTL in cell-mediated immune response.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.