



Recombinant Mouse Nidogen-1 (Nid1), partial

Product Code	CSB-BP015802MO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P10493
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	SPQRVNGKVKGRIFVGSSQVPVVFENTDLHSYVVMNHGRSYTAISTIPETVGY SLLPLAPIGGIIGWMFAVEQDGFKNQFSITGGEFTRQAEVTF LGHPGKLV LKQQ FSGIDEHGHLTISTELEGRVPQIPY GASVHIEPYTELYHYSSSVITSSSTREYTV MEPDQDGAAPSHTHIYQWRQTITFQECAHDDARPALPSTQQLSVDSVFLYN KEERILRYALSNSIGPVRDGPDA
Research Area	Others
Source	Baculovirus
Target Names	Nid1
Protein Names	Recommended name: Nidogen-1 Short name= NID-1 Alternative name(s): Entactin
Expression Region	428-665aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	This gene encodes a member of the nidogen family of basement membrane glycoproteins. The protein interacts with several other components of basement membranes, and may play a role in cell interactions with the extracellular matrix.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.