



Recombinant Chicken Collagen alpha-1 (XII) chain (COL12A1), partial

Product Code	CSB-EP005719CH-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P13944
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	DLVFLVDGWSVSGRNNFRYILDFMVALVSAFDIGEEKTRVGVVQYSSDTRTEF NLNQYFRRSDLLDAIKRIPYKGGNTMTGEAIDYLVKNTFTESAGARKGFPKVAI VITDGKAQDEVEIPARELRNIGVEVFSLGIKAADAKELKLIASQPSLKHVFNVAN FDGIVDIQNEI
Research Area	Others
Source	E.coli
Target Names	COL12A1
Protein Names	Recommended name: Collagen alpha-1(XII) chain Alternative name(s): Fibrochimerin
Expression Region	139-311aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	This gene encodes the alpha chain of type XII collagen, a member of the FACIT (fibril-associated collagens with interrupted triple helices) collagen family. Type XII collagen is a homotrimer found in association with type I collagen, an association that is thought to modify the interactions between collagen I fibrils and the surrounding matrix. Alternatively spliced transcript variants encoding different isoforms have been identified.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.