



Recombinant Bovine Pulmonary surfactant-associated protein B (SFTP B)

Product Code	CSB-EP021173BO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P15781
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	FPIPIPYCWL CRTLIKRIQAVIPKGV LAMTVAQVCHV VPLL VGGICQCLVERY SVI LLDTLLGRMLPQLVCGLVLRCS S
Research Area	Others
Source	E.coli
Target Names	SFTP B
Protein Names	Recommended name: Pulmonary surfactant-associated protein B Short name= SP-B Alternative name(s): 6 kDa protein Pulmonary surfactant-associated proteolipid SPL(Phe)
Expression Region	188-266aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.