




# Recombinant Mouse Serpin H1 (Serpinh1)

<b>Product Code</b>	CSB-EP021087MO
<b>Relevance</b>	Binds specifically to collagen. Could be involved as a chaperone in the biosynthetic pathway of collagen.
<b>Abbreviation</b>	Recombinant Mouse Serpinh1 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P19324
<b>Alias</b>	47 kDa heat shock protein Collagen-binding protein Short name: Colligin Serine protease inhibitor J6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	AEVKKPLEAAAPGTAEKLSKATTLAERSTGLAFSLYQAMAKDQAVENILLSPL VVASSLGLVSLGGKATTASQAKAVLSAEKLRDEEVHTGLGELLRSLSNSTARN VTWKLGSRLYGPSSVSFADDFVRSSKQHYNCEHSKINFRDKRSALQSINWA SQTTDGKLPEVTKDVERTDGALLVNAMFFKPHWDEKFHHKMVDNRGFMVTR SYTVGVMTMMHRTGLYNYYDDEKEKLMVEMPLAHLSSLIILMPHHVEPLERL EKLLTKEQLKAWMGKMQKKAVALSLPKGVVEVTHDLQKHLAAGLGLTEAIDKNK ADLSRMSGKKDLYLASVFHATAFEWDTEGNPFDQDIYGREELRSPKLFYADH PFIFLVRDNQSGSLLFIGRLVVRPKGDKMRDEL
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	Serpinh1
<b>Protein Names</b>	Recommended name: Serpin H1 Alternative name(s): 47 kDa heat shock protein Collagen-binding protein Short name= Colligin Serine protease inhibitor J6
<b>Expression Region</b>	18-417aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	60.8kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.