



Recombinant Human Nucleolin (NCL), partial

Product Code	CSB-MP015535HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P19338
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	VKLAKAGKNQGDPKKMAPPPKEVEEDSEDEEMSEDEEDDSSGEEVVIPQKKG KCAAATSAAKVVVSPTKKVAVATPAKKA AVTPGKKAATPAKKT VTPAKAVTT PGKKGATPGKALVATPGKKGAAIPAKGAKNGKNAKKEDSDEEEDDDSEEDDE DDEDEDEDEDEIEPAAMKAAAAAPASEDEDEDDEDDEDDEDDDDDEEDDSEEEA METTPAKGKKA AKVVPVAKNVAEDEDEEEDDEDDEDDEDDDEDDEDDEDDED EEEEEEEEEEPVKEAPGKRKKEMAKQKAAPEAKKQKVEGTEPTTAFNLFVGN LNFNKSAPELKTGISDVFAKNDLAVVDVRIGMTRKFGYVDFESAEDLEKALELT GLKVFGNEIKLEKPKGKDSKKERDARTLLAKNLPYKVTQDELKEVFEDAAEIRL VSKDGKSKGIAYIEFKTEADA EKTFEKQGT EIDGRSISLYYTGEK GQNQDYRG GKNSTWS
Research Area	Epigenetics and Nuclear Signaling
Source	Mammalian cell
Target Names	NCL
Protein Names	Recommended name: Nucleolin Alternative name(s): Protein C23
Expression Region	2-482aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.