



Recombinant Bovine Apolipoprotein A-I (APOA1), partial

Product Code	CSB-MP001913BO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P15497
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	DDPQSSWDRVKDFATVYVEAIKDSGRDYVAQFEASALGKQLNLKLLDNWDTL ASTLSKVREQLGPVTQEFWDNLEKETASLRQEMHKDLEEVKQKVPYLDEFQ KKWHEEVEIYRQKVAPLGEEFREGARQKVQELQDKLSPLAQELRDRARAHVE TLRQQLAPYSDDLRLRLTARLEALKEGGGSLAEYHAKASEQLKALGEKAKPVL EDLRQGLLPVLESLKVSILAAIDEASKKLNAAQ
Research Area	Others
Source	Mammalian cell
Target Names	APOA1
Protein Names	Recommended name: Apolipoprotein A-I Short name= Apo-AI Short name= ApoA-IA Alternative name(s): Apolipoprotein A1 Cleaved into the following chain: 1. Truncated apolipoprotein A-I
Expression Region	25-265aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	This gene encodes apolipoprotein A-I, which is the major protein component of high density lipoprotein (HDL) in plasma. The protein promotes cholesterol efflux from tissues to the liver for excretion, and it is a cofactor for lecithin cholesterolacyltransferase (LCAT) which is responsible for the formation of most plasma cholesteryl esters. This gene is closely linked with two other apolipoprotein genes on chromosome 11. Defects in this gene are associated with HDL deficiencies, including Tangier disease, and with systemic non-neuropathic amyloidosis.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.