



Recombinant Mouse Interleukin-2 receptor subunit beta (Il2rb), partial

Product Code	CSB-BP011650MO1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P16297
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	ASAAVKNCSHLECFYNSRANVSCMWSHEEALNVTTCVHAKSNLRHWNKTC ELTLVRQASWACNLILGSFPESQSLTSVDLLDINVCWEEKGWRRVKTCDFHP FDNLRLVAPHSLQVLHIDTQRCNISWKVSQVSHYIEPYLEFEARRRLLGHSWE DASVLSLKQRQQWLFLEMLIPSTSYEVQVRVKAQRNNTGTWSPWSQPLTFRT RPADPMKE
Research Area	Others
Source	Baculovirus
Target Names	Il2rb
Protein Names	Recommended name: Interleukin-2 receptor subunit beta Short name= IL-2 receptor subunit beta Short name= IL-2R subunit beta Short name= IL-2RBA Alternative name(s): High affinity IL-2 receptor subunit beta p70-75 CD_antigen=
Expression Region	24-240aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.