



# Recombinant Mouse Microtubule-associated protein 2 (Map2), partial

<b>Product Code</b>	CSB-MP013408MO
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P20357
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	GMEGEKLPPVPFAQTFGTNLEDRKQSTEPSIVMPSIGLSAEPPAPKEPKDWFIEMPTESSKDEWGLAAPISPGPLTPMREKDVLEDIPRWEGKQFDSMPMPSPFHGGSFTLPLDTMKNERVSEGPRPFAPVFFQSDDKVSLQDPSALATSKESSKDEEPLKDKADKVADVSISEVTTLLGNVHSPVVEGYVGENISGEVKVTTDQEKKETSAPSVQEPTLTETEPQTKLDEKSTVSIIEEAVAKKEESFKLRDDKTGVIQTSTEQSFSKEDQKGQEHTIDELKQDSFPISLEQAVTDAAMTSKTLGKVTSEPEAVSERREIQGLFEEKTADKNKLEGAGSATIAEVEMP
<b>Research Area</b>	Others
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Map2
<b>Protein Names</b>	Recommended name: Microtubule-associated protein 2 Short name= MAP-2
<b>Expression Region</b>	213-559aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Target Details</b>	This gene encodes a protein that belongs to the microtubule-associated protein family. The proteins of this family are thought to be involved in microtubule assembly, which is an essential step in neurogenesis. The products of similar genes in rat and mouse are neuron-specific cytoskeletal proteins that are enriched in dendrites, implicating a role in determining and stabilizing dendritic shape during neuron development. A number of alternatively spliced variants encoding distinct isoforms have been described.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.