



Recombinant Mouse Integrin alpha-L (Itgal), partial

Product Code	CSB-EP011875MO-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P24063
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	DLVFLFDGSQSLDRKDFEKILEFMKDVMRKLSNTSYQFAAVQFSTDCRTEFTF LDYVKQNKPNPDVLLGVSQPMFLLTNTFRAINYVVAHVFKESGARPDATKVLVI ITDGEASDKGNISAAHDITRYIIGIGKHFVSVQKQKTLHIFASEPVEEFVKILDTFE KDKDLFTDL
Research Area	Others
Source	E.coli
Target Names	Itgal
Protein Names	Recommended name: Integrin alpha-L Alternative name(s): CD11 antigen-like family member A Leukocyte adhesion glycoprotein LFA-1 alpha chain Short name= LFA-1A Leukocyte function-associated molecule 1 alpha chain Lymphocyte antigen
Expression Region	153-325aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	ITGAL encodes the integrin alpha L chain. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. This I-domain containing alpha integrin combines with the beta 2 chain (ITGB2) to form the integrin lymphocyte function-associated antigen-1 (LFA-1), which is expressed on all leukocytes. LFA-1 plays a central role in leukocyte intercellular adhesion through interactions with its ligands, ICAMs 1-3 (intercellular adhesion molecules 1 through 3), and also functions in lymphocyte costimulatory signaling. Two transcript variants encoding different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.