




Recombinant Mouse Protein Wnt-3a (Wnt3a)

Product Code	CSB-EP026136MO
Relevance	Ligand for members of the frizzled family of seven transmembrane receptors. Wnt-3 and Wnt-3a play distinct roles in cell-cell signaling during morphogenesis of the developing neural tube.
Abbreviation	Recombinant Mouse Wnt3a protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P27467
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	SYPIWWSLAVGPPQYSSLSTQPILCASIPGLVPKQLRFCRNYVEIMPSVAEGVKA GIQECQHQRGRRWNCCTTVSNSLAIFGPVLDKATRESAFVHAIASAGVAFVAVT RSCAEGSAAICGCSSRLQGSPGEGWKWGGCSEDIFFGGMVSREFADAREN PDARSAMNRHNNEAGRQAIASHMHLKCKCHGLSGSCEVKTCWWSQPDFRTI GDFLKD KYDSASEMVVEKHRESRGWVETLRPRYTYFKVPTERDLVYYEASPN FCEPNPETGSFGTRDRTC NVSSHGIDGCDLLCCGRGHNARTERRREKCHCV FHWCCYVSCQECTRVYDVHTCK
Research Area	Developmental Biology
Source	E.coli
Target Names	Wnt3a
Protein Names	Wnt-3a
Expression Region	19-352aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	40.9 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.